



# ETHANOL BREATH STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION  
Toxicology Section/Breath Alcohol Program  
SFN 59282 (5-2011)

Chemical Test Operator Name (Print) <i>Johnny Z. Lawson</i>		
Location <i>Bowman Co SD</i>	Intoxilyzer® Serial Number <i>80-004188</i>	
Gas Lot Number <i>B/6512080A2</i>	Gas Cylinder Number <i>31</i>	Gas Expiration Date <i>8/01/2014</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).

Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).

- A. 0. 079 AC
- B. 0. 078 AC
- C. 0. 079 AC

- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Breath Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>Johnny Z. Lawson</i>	Date <i>12-13-13</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Kasher</i>	Date <i>03 Jan 14</i>

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004188  
Location = BOWM      8164.13.00 06/09  
12/13/2013      11:29

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:30
02 Std. Gas	0.079	11:30
03 Room Air	0.000	11:30
04 Std. Gas	0.078	11:31
05 Room Air	0.000	11:31
06 Std. Gas	0.079	11:32
07 Room Air	0.000	11:32


Lot No = 16512080A2

Cyl No = 31

Exp Date = 08/01/2014

County = 06

Oper No. = 020814

  
Operator Signature  
JOHNNY LAWSON

Remarks:

Form 106-I8000