



ETHANOL BREATH STANDARD CYLINDER REPORT
OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (1-2010)

Only Chemical Test Operators may change the Ethanol Breath Standard and complete this form.

Chemical Test Operator Name: <u>David Feland</u>
Location: <u>Kidder So.</u>
Intoxilyzer® Serial Number: <u>80-004186</u>
Gas Lot Number: <u>659357</u> Gas Cylinder Number: <u>37</u> Expiration Date: <u>09/03/11</u>
<input type="checkbox"/> New Cylinder Installed <input checked="" type="checkbox"/> Cylinder Not Changed
Instrument Test Results (Report to 3 Digits; Example, 0.081%) Test 1: <u>0.082</u> % Test 2: <u>0.082</u> % Test 3: <u>0.082</u> % <p style="text-align: right;">(ATTACH TEST RECORD)</p>

David Feland
Chemical Test Operator Signature

16 Dec 10
Date

Deb Kashner
Reviewed By

22 Dec 2010
Date

80-004186
p. 95 of 96
AM

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004186
Location = STEE 8164.13.00 06/09
12/16/2010 14:57

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:57
02 Std. Gas	0.082	14:58
03 Room Air	0.000	14:58
04 Std. Gas	0.082	14:59
05 Room Air	0.000	14:59
06 Std. Gas	0.082	14:59
07 Room Air	0.000	15:00

Lot No = 659357
Cyl No = 37
Exp Date = 09/03/2011
County = 22 Oper No. = 119404



Operator Signature
DAVID FELAND

Remarks:

Form 106-I8000

80-004186
P. 96 of 96
FM