



ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program

SFN 59282 (10-2013)

Chemical Test Operator Name (Print) THOMAS HERZIG		
Location STURMAN LEC (JAIL)	Intoxilyzer® Serial Number 80-005-361	
Gas Lot Number 19615080AS	Gas Cylinder Number 034	Gas Expiration Date 09/05/2010

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0. 081 AC
 - B. 0. 081 AC
 - C. 0. 080 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature 	Date 12/31/15
Reviewed By (Crime Laboratory Use Only) Deb Shanaver	Date 11 Jan 16 2016

*JAIL UNIT
New CYLINDER
FRONTAL*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005361
Location = JAME 8164.13.00 06/09
12/31/2015 16:46

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:48
02 Std. Gas	0.081	16:48
03 Room Air	0.000	16:48
04 Std. Gas	0.081	16:49
05 Room Air	0.000	16:49
06 Std. Gas	0.080	16:50
07 Room Air	0.000	16:50


Lot No = 19615080A5

Cyl No = 34

Exp Date = 09/05/2017

County = 47

Oper No. = 039706



Operator Signature
THOMAS HERZIG

Remarks:

Form 106-I8000

JAIL UNIT
NEW CYLINDER
INSTALL