

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-005940

Location: Griggs Co SO

Check When Done

- A. Is the breath tube heated? ✓
- B. Is the diagnostic check complete? ✓
- C. Is the time, date, and year correct (re-set if necessary)? ✓
- D. Print test. (Attach test record.) ✓
- E. Low AC. Use  $\leq 0.03$  AC in ACA mode. (Attach test record.)  
Sim SN DR5113 Lot No. 201507F Rep AC 0.020 ✓
- F. High AC. Use  $\geq 0.25$  AC in ACA mode. (Attach test record.)  
Sim SN DR3378 Lot No. 16060 Rep AC 0.300 ✓
- G. Interferent check. Use 0.05% acetone in ethanol solution in ABA mode.  
(Attach test record.)  
Sim SN DR3803 Lot No. ILS4 Solution 0.10 AC + 0.05% acetone ✓
- H. RFI check. Run CMS mode. Key radio during CMS test.  
(Attach test record.) ✓
- I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.  
(Attach test record.) 0.080 AC

Lot No. 19615080A5 Cyl No. 19 Exp. Date 09/05/2017  
Test 1 0.079 Test 4 0.078 Test 7 0.079  
Test 2 0.078 Test 5 0.078 Test 8 0.079  
Test 3 0.079 Test 6 0.079 Test 9 0.079  
Average 0.079

J. Tank pressure, Level 3 Func D Sub G P: 246 psi Reg 250 psi

K. Remarks/Maintenance record: Annual inspection

09 May 17  
Date

10 MAY 17  
Date

Deb Shanover  
Field Inspector's Signature

Charles E. Egan  
Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005940  
Location = COOP                8164.13.00 06/09  
05/09/2017                      15:17

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?         1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            19615080A5  
Standard Cyl #?            19  
Standard Expiration?       09/05/2017  
Oper No?                    777777

Flow Cal. Date:            06/21/2013  
Slope                      688  
Intercept                  -659053

IR Calibration Date:        04/02/2015  
                                    3um                      9um

	3um	9um
0th Coef(*100):	-18917	-25088
1st Coef(*100):	267616	136158
2nd Coef(*100):	3374	1493
H2O adj(mg/l*10k):	427	417

\*\*\*\*\* Printer Test End \*\*\*\*\*

Deb Shanaver

Operator Signature  
DEB SHANAVER

Remarks:

Print test

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005940  
Location = COOP                 8164.13.00 06/09  
05/09/2017                        15:32

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:33
02 Std. Sol.	0.019	15:33
03 Room Air	0.000	15:34
04 Std. Sol.	0.019	15:35
05 Room Air	0.000	15:35
06 Std. Sol.	0.019	15:36
07 Room Air	0.000	15:36

08 Sim Temp = 34.0°C

Simul Ser No = DR5113

Std Sol No = 201507F

County = 20

Oper No. = 777777

Deb Shanaver

Operator Signature  
DEB SHANAVER

Remarks:

Low AC 0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005940  
Location = COOP      8164.13.00 06/09  
05/09/2017      15:38

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:39
02 Std. Sol.	0.312	15:39
03 Room Air	0.000	15:40
04 Std. Sol.	0.310	15:41
05 Room Air	0.000	15:41
06 Std. Sol.	0.309	15:42
07 Room Air	0.000	15:42

08 Sim Temp = 34.0°C

Simul Ser No = DR3378

Std Sol No = 16060

County = 20

Oper No. = 777777



Operator Signature  
DEB SHANAVER

Remarks:

*High AC 0.300 AC*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005940  
Location = COOP                 8164.13.00 06/09  
05/09/2017                        15:46

Test	AC	Time
01 Room Air	0.000	15:47
02 *Subject Test	INT*	15:47
03 Room Air	0.000	15:48

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                 Weight = 150  
Test = DUI                         Cit = NA  
Dr. Lic. = ND/DIS921456  
Lot No = 19615080A5  
Cyl No = 19  
Expiration Date = 09/05/2017  
County = 20                        Oper No. = 777777

Deb Shanaver  
Operator Signature  
DEB SHANAVER

Remarks:

*Interferent detected check*

Form 106-I8000    0.10A ct + 0.05%  
acetone

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005940  
Location = COOP                8164.13.00 06/09  
05/09/2017                      15:30

Test	AC	Time
01 Diagnostic	OK	15:31
02 Room Air	RFI*	15:31
03 Room Air	0:000	15:31

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                Weight = 150  
Test = DUI                        Cit = NA  
Dr. Lic. = ND/DIS921456  
Lot No = 19615080A5  
Cyl No = 19  
Expiration Date = 09/05/2017  
County = 20                      Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Deb Shanaver  
Operator Signature  
DEB SHANAVER

Remarks:

RFI check

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005940  
Location = COOP      8164.13.00 06/09  
05/09/2017      15:18

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:19
02 Std. Gas	0.079	15:19
03 Room Air	0.000	15:19
04 Std. Gas	0.078	15:20
05 Room Air	0.000	15:20
06 Std. Gas	0.079	15:21
07 Room Air	0.000	15:21

Lot No = 19615080A5  
Cyl No = 19  
Exp Date = 09/05/2017  
County = 20      Oper No. = 777777

Deb Shanaver  
Operator Signature  
DEB SHANAVER

Remarks:

Calibration check 0.080 AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-005940  
Location = COOP                    8164.13.00 06/09  
05/09/2017                            15:22

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:22
02 Std. Gas	0.078	15:23
03 Room Air	0.000	15:23
04 Std. Gas	0.078	15:23
05 Room Air	0.000	15:24
06 Std. Gas	0.079	15:24
07 Room Air	0.000	15:25

Lot No = 19615080A5  
Cyl No = 19  
Exp Date = 09/05/2017  
County = 20                            Oper No. = 777777

Deb Shanaver  
Operator Signature  
DEB SHANAVER

Remarks:

Calibration check 0.080AC

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-005940  
Location = COOP                    8164.13.00 06/09  
05/09/2017                            15:25

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:26
02 Std. Gas	0.079	15:26
03 Room Air	0.000	15:27
04 Std. Gas	0.079	15:27
05 Room Air	0.000	15:27
06 Std. Gas	0.079	15:28
07 Room Air	0.000	15:28

Lot No = 19615080A5  
Cyl No = 19  
Exp Date = 09/05/2017  
County = 20                            Oper No. = 777777

  
\_\_\_\_\_  
Operator Signature  
DEB SHANAVER

Remarks:

*Calibration check 0.080AC*