

INTOXILYZER® 8000 INSPECTION

Serial No.: 80-006495

Location: Fargo PD - Dept. of Health Bldg

- Check When Done
- A. Is the breath tube heated? ✓
 - B. Is the diagnostic check complete? ✓
 - C. Is the time, date, and year correct (re-set if necessary)? ✓
 - D. Print test. (Attach test record.) ✓
 - E. Low AC. Use ≤ 0.03 AC in ACA mode. (Attach test record.)
Sim SN DR5113 Lot No. 201507F Rep AC 0.020 ✓
 - F. High AC. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN DR 3378 Lot No. 16060 Rep AC 0.300 ✓
 - G. Interferent check. Use 0.05% acetone in ethanol solution in ABA mode.
(Attach test record.)
Sim SN DR 3803 Lot No. ICS 4 Solution 0.10AC + 0.05% acetone ✓
 - H. RFI check. Run CMS mode. Key radio during CMS test.
(Attach test record.) ✓
 - I. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.) 0.080 AC

Lot No. 09216080A3 Cyl No. 37 Exp. Date 05/05/2018

Test 1 <u>0.082</u>	Test 4 <u>0.083</u>	Test 7 <u>0.082</u>
Test 2 <u>0.082</u>	Test 5 <u>0.082</u>	Test 8 <u>0.082</u>
Test 3 <u>0.083</u>	Test 6 <u>0.083</u>	Test 9 <u>0.083</u>

Average 0.082

J. Tank pressure, Level 3 Func D Sub G P: 842 psi Reg 860 psi

K. Remarks/Maintenance record: Annual inspection

09 May 17
Date

10 MAY 17
Date

Deb Shanaver
Field Inspector's Signature

Charles E. Edin
Reviewed by

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006495
Location = FAPD 8164.13.00 06/09
05/09/2017 09:26

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 09216080A3
Standard Cyl #? 37
Standard Expiration? 05/05/2018
Oper No? 777777

Flow Cal. Date: 08/18/2015
Slope 678
Intercept -657244

IR Calibration Date: 08/18/2015
 3um 9um

0th Coef(*100): -22405 -18664
1st Coef(*100): 270469 136654
2nd Coef(*100): 3040 1360
H2O adj(mg/l*10k): 699 443

***** Printer Test End *****

Deb Shanaver

Operator Signature
DEB SHANAVER

Remarks:

Print test

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006495
Location = FAPD 8164.13.00 06/09
05/09/2017 09:39

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:40
02 Std. Sol.	0.020	09:41
03 Room Air	0.000	09:41
04 Std. Sol.	0.020	09:42
05 Room Air	0.000	09:43
06 Std. Sol.	0.020	09:43
07 Room Air	0.000	09:44

08 Sim Temp = 34.0°C

Simul Ser No = DR5113
Std Sol No = 201507F
County = 09 Oper No. = 777777



Operator Signature
DEB SHANAVER

Remarks:

Low AC 0.020AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006495
Location = FAPD 8164.13.00 06/09
05/09/2017 09:45

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:45
02 Std. Sol.	0.310	09:46
03 Room Air	0.000	09:47
04 Std. Sol.	0.309	09:48
05 Room Air	0.000	09:48
06 Std. Sol.	0.309	09:49
07 Room Air	0.000	09:50

08 Sim Temp = 34.0°C

Simul Ser No = DR3378

Std Sol No = 16060

County = 09

Oper No. = 777777

Deb Shanaver

Operator Signature
DEB SHANAVER

Remarks:

High AC 0.300 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006495
Location = FAPD 8164.13.00 06/09
05/09/2017 09:55

Test	AC	Time
01 Room Air	0.000	09:55
02 *Subject Test	INT*	09:56
03 Room Air	0.000	09:57

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 09216080A3
Cyl No = 37
Expiration Date = 05/05/2018
County = 09 Oper No. = 777777

Deb Shanaver

Operator Signature
DEB SHANAVER

Remarks:

Interferent detected check

Form 106-I8000

*0.10AC + 0.05^gro
acetone*

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006495
Location = FAPD 8164.13.00 06/09
05/09/2017 09:58

Test	AC	Time
01 Diagnostic	OK	09:59
02 Room Air	RFI*	09:59
03 Room Air	0.000	09:59

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 09216080A3
Cyl No = 37
Expiration Date = 05/05/2018
County = 09 Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Deb Shanaver

Operator Signature
DEB SHANAVER

Remarks:

RFI check

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006495
Location = FAPD 8164.13.00 06/09
05/09/2017 09:27

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:28
02 Std. Gas	0.082	09:28
03 Room Air	0.000	09:29
04 Std. Gas	0.082	09:29
05 Room Air	0.000	09:30
06 Std. Gas	0.083	09:30
07 Room Air	0.000	09:31

Lot No = 09216080A3

Cyl No = 37

Exp Date = 05/05/2018

County = 09

Oper No. = 777777

Deb Shanaver

Operator Signature
DEB SHANAVER

Remarks:

Calibration check 0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006495
Location = FAPD 8164.13.00 06/09
05/09/2017 09:31

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:32
02 Std. Gas	0.083	09:32
03 Room Air	0.000	09:33
04 Std. Gas	0.082	09:33
05 Room Air	0.000	09:34
06 Std. Gas	0.083	09:34
07 Room Air	0.000	09:35

Lot No = 09216080A3
Cyl No = 37
Exp Date = 05/05/2018
County = 09 Oper No. = 777777

Deb Shanaver

Operator Signature
DEB SHANAVER

Remarks:

Calibration check 0.080 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006495
Location = FAPD 8164.13.00 06/09
05/09/2017 09:35

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:35
02 Std. Gas	0.082	09:36
03 Room Air	0.000	09:36
04 Std. Gas	0.082	09:37
05 Room Air	0.000	09:37
06 Std. Gas	0.083	09:38
07 Room Air	0.000	09:38

Lot No = 09216080A3
Cyl No = 37
Exp Date = 05/05/2018
County = 09 Oper No. = 777777

Deb Shanaver

Operator Signature
DEB SHANAVER

Remarks:

Calibration check 0.080 AC

Form 106-I8000