



**ETHANOL GAS STANDARD CYLINDER REPORT**  
 OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION  
 Toxicology Section/Breath Alcohol Program  
 SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Tracy L. Meidinger</i>		
Location <b>UNDP</b>	Intoxilyzer® Serial Number <b>80-005357</b>	
Gas Lot Number <b>03417080A3</b>	Gas Cylinder Number <b>32</b>	Gas Expiration Date <b>3-05-2019</b>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- REN* 2. Perform an ACA Test (Level 1, Function C).  
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
  - A. 0.081 AC
  - B. 0.082 AC
  - C. 0.083 AC *REN*
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>Tracy L. Meidinger</i>	Date <b>9-11-2017</b>
Reviewed By (Crime Laboratory Use Only) <i>Roberta Luque-Nemmo</i>	Date <b>9-13-2017</b>


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-005357  
Location = UNDP                    8164.14.00 09/16  
09/11/2017                            06:39

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	06:39
02 Std. Gas	0.081	06:40
03 Room Air	0.000	06:40
04 Std. Gas	0.082	06:41
05 Room Air	0.000	06:41
06 Std. Gas	0.082	06:42
07 Room Air	0.000	06:42

Lot No = 03417080A3  
Cyl No = 32  
Exp Date = 03/05/2019  
County = 18                            Oper No. = 119409

  
\_\_\_\_\_  
Operator Signature  
TRACY MEIDINGER

Remarks: *Cylinder Install #32*