



**ETHANOL GAS STANDARD CYLINDER REPORT**  
 OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION  
 Toxicology Section/Breath Alcohol Program  
 SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Alex Schwab</i>		
Location <i>Powers Lake Police Dept.</i>	Intoxilyzer® Serial Number <i>88-005942</i>	
Gas Lot Number <i>03417080A3</i>	Gas Cylinder Number <i>5</i>	Gas Expiration Date <i>03/05/2019</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
  - A. 0.080 AC
  - B. 0.079 AC
  - C. 0.080 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>Alex Schwab</i>	Date <i>9-22-17</i>
Reviewed By (Crime Laboratory Use Only) <i>Roberta Suega-Nunne</i>	Date <i>9-22-17</i>

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005942  
Location = PLPD      8164.13.00 06/09  
09/22/2017      11:43

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:44
02 Std. Gas	0.080	11:44
03 Room Air	0.000	11:45
04 Std. Gas	0.079	11:45
05 Room Air	0.000	11:46
06 Std. Gas	0.080	11:46
07 Room Air	0.000	11:47

Lot No = 03417080A3

Cyl No = 5

Exp Date = 03/05/2019

County = 07

Oper No. = 131308



Operator Signature  
ALEX SCHWAB

Remarks:

Form 106-I8000