



# ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION  
Toxicology Section/Breath Alcohol Program  
SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Gabriel Johnston</i>		
Location <i>Cavalier County SO</i>	Intoxilyzer® Serial Number <i>80-004938</i>	
Gas Lot Number <i>29616080A2</i>	Gas Cylinder Number <i>46</i>	Gas Expiration Date <i>09-05-2018</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
  - A. 0.082 AC
  - B. 0.082 AC
  - C. 0.082 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>[Signature]</i>	Date <i>10-26-2017</i>
Reviewed By (Crime Laboratory Use Only) <i>Roberto Aguero-Munoz</i>	Date <i>10-30-2017</i>


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004938  
Location = LANG      8164.14.00 09/16  
10/26/2017      13:40

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:41
02 Std. Gas	0.082	13:41
03 Room Air	0.000	13:42
04 Std. Gas	0.082	13:42
05 Room Air	0.000	13:42
06 Std. Gas	0.082	13:43
07 Room Air	0.000	13:43

Lot No = 20016080A2  
Cyl No = 46  
Exp Date = 09/05/2018  
County = 10      Oper No. = 021009

  
\_\_\_\_\_  
Operator Signature  
GABRIEL JOHNSTON

Remarks: *ACA TEST*