



ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Nathan C Moen</i>		
Location <i>Grand Forks County Correctional Center</i> ^{BFSO}		Intoxilyzer® Serial Number <i>80-006501</i>
Gas Lot Number <i>19817080A4</i>	Gas Cylinder Number <i>17</i>	Gas Expiration Date <i>09/05/2019</i>

Check When Done:

1. Scan/Enter Gas Cylinder Information (Level 1, Function S).

2. Perform an ACA Test (Level 1, Function C).

*ACA test was done.
BFSO
Crime Lab*

Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).

- A. *0.081* AC
- B. *0.082* AC
- C. *0.084* AC

3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.

4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.

5. Send the Following to the Crime Laboratory:

- A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
- B. ACA Test Record.

Chemical Test Operator Signature <i>Nathan C Moen</i>	Date <i>11/22/2017</i>
Reviewed By (Crime Laboratory Use Only) <i>Roberta Hugu-Nimmo</i>	Date <i>11/30/2017</i>


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006501
Location = GFSO 8164.14.00 09/16
11/22/2017 21:33

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	21:33
02 Std. Gas	0.081	21:34
03 Room Air	0.000	21:34
04 Std. Gas	0.082	21:35
05 Room Air	0.000	21:35
06 Std. Gas	0.081	21:36
07 Room Air	0.000	21:36

Lot No = 19817080A4
Cyl No = 17
Exp Date = 09/05/2019
County = 18 Oper No. = 100452


Operator Signature
NATHAN MOEN

Remarks: