



ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (10-2013)

Chemical Test Operator Name (Print) Daniel Garvens		
Location CAVA		Intoxilyzer® Serial Number 80-006491
Gas Lot Number 19817080A4	Gas Cylinder Number 3	Gas Expiration Date 9-5-19

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0. 0 8 2 AC
 - B. 0. 0 8 2 AC
 - C. 0. 0 8 2 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>[Signature]</i> 4502	Date 1-4-18
Reviewed By (Crime Laboratory Use Only) <i>[Signature]</i>	Date 1-4-18

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006491
Location = CAVA 8164.14.00 09/16
01/04/2018 10:51

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:51
02 Std. Gas	0.082	10:52
03 Room Air	0.000	10:52
04 Std. Gas	0.082	10:53
05 Room Air	0.000	10:53
06 Std. Gas	0.082	10:54
07 Room Air	0.000	10:54

Lot No = 19817080A4
Cyl No = 3
Exp Date = 09/05/2019
County = 34 Oper No. = 060849



Operator Signature
DANIEL GARVENS

Remarks: *New Gas Cylinder Install*