



ETHANOL GAS STANDARD CYLINDER REPORT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59282 (10-2013)

Chemical Test Operator Name (Print) FRANK EBERLE		
Location BOWMAN CO. Sheriff's Office	Intoxilyzer® Serial Number 80-006496	
Gas Lot Number 19817080A4	Gas Cylinder Number 5	Gas Expiration Date 09-05-2019

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0.082 AC
 - B. 0.082 AC
 - C. 0.083 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature Frank Eberle	Date 1-22-2018
Reviewed By (Crime Laboratory Use Only) Roberta August Nimmo	Date 1-24-2018


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006496
Location = BOWM 8164.13.00 06/09
01/22/2018 11:38

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:39
02 Std. Gas	0.082	11:39
03 Room Air	0.000	11:40
04 Std. Gas	0.082	11:40
05 Room Air	0.000	11:41
06 Std. Gas	0.083	11:41
07 Room Air	0.000	11:42

Lot No = 19817080A4
Cyl No = 5
Exp Date = 09/05/2019
County = 06 Oper No. = 130015



Operator Signature
FRANK EBERLE

Remarks:

Form 106-I8000