



ETHANOL GAS STANDARD CYLINDER REPORT
OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (10-2013)

Chemical Test Operator Name (Print) KEVIN ELDRIDGE		
Location DICKINSON LEC	Intoxilyzer® Serial Number 80-005947	
Gas Lot Number 19817080A4	Gas Cylinder Number 10	Gas Expiration Date 09/05/2019

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).

Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).

- A. 0.079 AC
- B. 0.079 AC
- C. 0.079 AC

- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature	Date 4/3/18
Reviewed By (Crime Laboratory Use Only)	Date 4/4/18


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005947
Location = DKPD 8164.14.00 09/16
04/03/2018 11:25

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:26
02 Std. Gas	0.079	11:26
03 Room Air	0.000	11:27
04 Std. Gas	0.079	11:27
05 Room Air	0.000	11:28
06 Std. Gas	0.079	11:28
07 Room Air	0.000	11:29

Lot No = 19817080A4
Cyl No = 10
Exp Date = 09/05/2019
County = 45 Oper No. = 130935


Operator Signature
KEVIN ELDRIDGE

Remarks: