



**ETHANOL GAS STANDARD CYLINDER REPORT**  
 OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION  
 Toxicology Section/Breath Alcohol Program  
 SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Darcy W Aberle</i>		
Location <i>Williams Co LEC</i>	Intoxilyzer® Serial Number <i>80-005946</i>	
Gas Lot Number <i>34917080A3</i>	Gas Cylinder Number <i>17</i>	Gas Expiration Date <i>02-05-2020</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).  
 A. *0.079* AC  
 B. *0.078* AC  
 C. *0.078* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>Darcy Aberle</i>	Date <i>07/15/2018</i>
Reviewed By (Crime Laboratory Use Only) <i>Roberta Ruggen-Nimmo</i>	Date <i>01/22/2019</i>


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005946  
Location = WILL      8164.14.00 09/16  
07/15/2018      16:23

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:24
02 Std. Gas	0.079	16:24
03 Room Air	0.000	16:25
04 Std. Gas	0.078	16:25
05 Room Air	0.000	16:26
06 Std. Gas	0.078	16:26
07 Room Air	0.000	16:26

Lot No = 34917080A3  
Cyl No = 17  
Exp Date = 02/05/2020  
County = 53      Oper No. = 010101

  
Operator Signature  
DARCY ABERLE

Remarks:

Form 106-I8000