



ETHANOL GAS STANDARD CYLINDER REPORT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Melby, Robert M.</i>		
Location <i>"CROS" Divide Co. SD.</i>	Intoxilyzer® Serial Number <i>80-005957</i>	
Gas Lot Number <i>34917080A3</i>	Gas Cylinder Number <i>32</i>	Gas Expiration Date <i>02-05-2020</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0.*078* AC
 - B. 0.*078* AC
 - C. 0.*078* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Robert M. Melby</i>	Date <i>08-23-2018</i>
Reviewed By (Crime Laboratory Use Only) <i>Robert M. Melby</i>	Date <i>08-24-2018</i>


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005957
Location = CROS 8164.14.00 09/16
08/23/2018 13:36

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:38
02 Std. Gas	0.078	13:38
03 Room Air	0.000	13:39
04 Std. Gas	0.078	13:39
05 Room Air	0.000	13:39
06 Std. Gas	0.078	13:40
07 Room Air	0.000	13:40

Lot No = 34917080A3
Cyl No = 32
Exp Date = 02/05/2020
County = 12 Oper No. = 100011




Operator Signature
ROBERT MELBY

Remarks:

Test Good

Form 106-I8000

The undersigned, having custody of the original record, certifies that the information hereon is a true and correct copy of the original document maintained as part of the activities of this office.



Signature

8-23-2018

Date