



ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program

SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Gabriel Johnston</i>		
Location <i>Cavalier County - Langdon</i>	Intoxilyzer® Serial Number <i>80-004938</i>	
Gas Lot Number <i>34917080A3</i>	Gas Cylinder Number <i>30</i>	Gas Expiration Date <i>02-05-2020</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0. 080 AC
 - B. 0. 080 AC
 - C. 0. 080 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>[Signature]</i>	Date <i>04-01-2018</i>
Reviewed By (Crime Laboratory Use Only) <i>[Signature]</i>	Date <i>09-05-2018</i>

Intoxilyzer Test Record and Checklist
NDDAG Crime Lab. Div., Bismarck, ND 58501

CHI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004938
Location = LANG 8164.14.00 09/16
09/01/2018 21:10

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	21:11
02 Std. Gas	0.080	21:11
03 Room Air	0.000	21:12
04 Std. Gas	0.080	21:12
05 Room Air	0.000	21:13
06 Std. Gas	0.080	21:13
07 Room Air	0.000	21:14

Lot No = 34917080A3
Cyl No = 30
Exp Date = 02/05/2020
County = 10 Oper No. = 021009

Operator Signature
GABRIEL JOHNSTON

Remarks: *NEW GAS CYLINDER*
[Signature] #4228
Form 106-18000