OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

INTOXILYZER® 8000 INSPECTION

Intoxilyzer® 8000 Serial Number: 80-0	006489	Location: GRAND FORKS P.D. (GFPD)
A. General Setup and Checks: 1. X Diagnostic Tests Pass at 2. X Breath tube heated 3. X Date, time and location 4. X Print test (Level 1,P). Stank monitor (Level 3,D a. Display: 409 psi F	code (Level 2,E). R sign and attach test r	eady Mode le-set if necessary. record. egulator within 50 psi.
Sim. Ser #: DR337 b.	Iowing tests (Level 1 C): Known Value ≤ Lot #: 17160 C): Known Value ≥ Lot #: 17350 The following tests I 1,B): Known Value #:	0.03 AC: <u>0.026</u> AC _Exp. Date: <u>6.8.19</u> : 0.25 AC: <u>0.300</u> AC _Exp. Date: <u>10.11.19</u> (Level 1,S). :: <u>0.10</u> AC + 0.05% Acetone Date: <u>NA</u> ads "RFI Detect" /alue <u>0.080</u> AC _Exp. Date: <u>2.5.20</u> Test 7 <u>0.080</u> Test 8 <u>0.080</u>
C. Remarks/Maintenance: ANNI	IAL INSPEC	TON
Instrument is acceptable to be used in	the field.	
Breath Analyst Signature NA		
Reviewed by		Date

02/18

Issuing Authority - CEE

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006489
Location = GFPD 8164.14.00 09/16
05/13/2019 19:26

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 34917080A3 Standard Cyl #? 12 02/05/2020 Standard Expiration? Oper No? 666666

Flow Cal. Date: 08/18/2015 Slope 663 Intercept -631787

IR Calibration Date: 08/18/2015
3um 9um

Oth Coef(*100): -24626 -16551 1st Coef(*100): 269218 130789 2nd Coef(*100): 3390 1656 H2O adj(mq/l*10k): 742 460

****** Printer Test End ********

Operator Signature CHARLES EDER

Remarks: PRINT TEST

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = GFPD 05/13/2019 Alcohol Analyzer SN 80-006489 8164.14.00 09/16 19:39

WET CAL CHECK

Test	AC	Time	
01 Room Air	0.000	19:39	
02 Std. Sol.	0.019	19:40	
03 Room Air	0.000	19:40	
04 Std. Sol.	0.019	19:41	
05 Room Air	0.000	19:42	
06 Std. Sol.	0.019	19:43	
07 Room Air	0.000	19:43	

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = DR3375 Std Sol No = 17160 County = 18

Oper No. = 666666

Operator Signature CHARLES EDER

LOW AC

Remarks:

0-020AC

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = GFPD 05/13/2019 Alcohol Analyzer SN 80-006489 8164.14.00 09/16 19:44

WET CAL CHECK

Test		AC	Time	
01 Room	Air	0.000	19:45	
02 Std.	Sol.	0.303	19:45	
03 Room	Air	0.000	19:46	
04 Std.	Sol.	0.303	19:47	
05 Room	Air	0.000	19:47	
06 Std.	Sol.	0.303	19:48	
07 Room	Air	0.000	19:49	

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = DR7111 Std Sol No = 17350

County = 18

Oper No. = 666666

Operator Signature CHARLES EDER

HIGH AC

0.300

Remarks:

Alcohol Analyzer CMI, Inc. Intoxilyzer North Dakota Model 8000 SN 80-006489 Location = GFPD 8164.14.00 09/16 05/13/2019

Test	t AC	
01 Room Air	0.000	19:24
02 *Subject Test	INT*	19:24
03 Room Air	0.000	19:25

19:23

Interferent Detected

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 01/01/1982

Sub Sex = FemaleWeight = 150

Cit = INTERFERENT Test = DUI Dr. Lic. = ND/DIS821456

Lot No = 34917080A3

Cyl No = 12

Expiration Date = 02/05/2020

County = 18Oper No. = 666666

> Operator Signature CHARLES EDER

Remarks:

Interferent Chack 0.10AC + 0.05% ACETONE

^{*}Invalid Test

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006489
Location = GFPD 8164.14.00 09/16
05/13/2019 19:31

Test	AC	Time	
01 Diagnostic	OK	19:32	
02 Room Air	0.000	19:33	
03 *Subject Test	RFI*	19:33	
04 Room Air	0.000	19:33	

*Invalid Test Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 01/01/1982

Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK

Dr. Lic. = ND/DIS821456

Lot No = 34917080A3

Cyl No = 12

Expiration Date = 02/05/2020

County = 18 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = GFPD
05/13/2019

Alcohol Analyzer SN 80-006489 8164.14.00 09/16 19:26

DRY CAL CHECK

Test	AC	Time	
01 Room Air	0.000	19:27	
02 Std. Gas	0.079	19:27	
03 Room Air	0.000	19:28	
04 Std. Gas	0.080	19:28	
05 Room Air	0.000	19:29	
06 Std. Gas	0.080	19:29	
07 Room Air	0.000	19:30	

Lot No = 34917080A3

Cyl No = 12

Exp Date = 02/05/2020

County = 18

Oper No. = 666666

Operator Signature

Remarks: CACIBRATION CHECK

0.080AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = GFPD
05/13/2019

Alcohol Analyzer SN 80-006489 8164.14.00 09/16 19:34

DRY CAL CHECK

Test			AC		Time	
	01	Room	Air	0	.000	19:35
	02	Std.	Gas	0	.080	19:35
	03	Room	Air	0	.000	19:36
	04	Std.	Gas	0	.081	19:36
	05	Room	Air	0	.000	19:37
	06	Std.	Gas	0	.080	19:37
	07	Room	Air	0	.000	19:38

Lot No = 34917080A3

Cyl No = 12

Exp Date = 02/05/2020

County = 18

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks: CACIBRATION CHECK

0.080 AC

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = GFPD 05/13/2019

Alcohol Analyzer SN 80-006489 8164.14.00 09/16 19:50

DRY CAL CHECK

Test			AC	Time	
	01	Room	Air	0.000	19:51
	02	Std.	Gas	0.080	19:51
	03	Room	Air	0.000	19:52
	04	Std.	Gas	0.080	19:52
	05	Room	Air	0.000	19:53
	06	Std.	Gas	0.081	19:53
	07	Room	Air	0.000	19:54

Lot No = 34917080A3

Cyl No = 12

Exp Date = 02/05/2020

County = 18

Oper No. = 666666

Remarks:

Operator Signature
CHARLES EDER

ALIBRATION CHECK 0.090 AC