OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program Intoxilyzer® 8000 Inspection

INTOXILYZER® 8000 INSPECTION

ntoxilyzer® 8000 Serial Number: 80-004193 Location: ToxL
 General Setup and Checks: 1.
 Tests (Sign and attach test records):
Configure dry gas standard for the following tests (Level 1,S).
2. Interferent Check (Level 1,B): Known Value: 0-10 AC + 0.05% Acetone Sim. Ser #: DR3451 Lot #: FCS 5 Exp. Date: NA Display reads "Interferent Detect" 3. Interferent Check (CMS Mode) Interferent Detect" 4. Interferent Check (CMS Mode) Interferent Detect" 4. Interferent Check (Level 1,C): Exp. Date: National Detect (CMS Mode) Interferent Detect (CMS Mode) Inte
Remarks/Maintenance: ANNUAL FUSPECTION, DURING CALIBRATION AD TUST MENT THE AUTO RANGE VALUE = 3, THIS BELOW 4. THIS INSTRUMENT WILL BE USED AS CLASSROOM ONLY UNTIL ANALYTICAL BENCH CAN BE UPDATED.
Instrument is acceptable to be used in the field.
Cambre Estal 6.19.2019
Breath-Analyst Signature NA Date
Reviewed by Date

02/18 Issuing Authority - CEE

SN 80-004193

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004193
Location = TOXL 8164.14.00 09/16
06/18/2019 16:48

******** Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer (Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 34418080A2 Standard Cyl #? 40 02/05/2021 Standard Expiration? Oper No? 666666

Flow Cal. Date: 06/18/2019 Slope 622 Intercept -723749

IR Calibration Date: 06/18/2019
3um 9um

Oth Coef(*100): -14782 -15525
1st Coef(*100): 267622 134211
2nd Coef(*100): 2469 1112
H2O adj(mg/l*10k): 264 225

******* Printer Test End ********

Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
06/18/2019

Alcohol Analyzer SN 80-004193 8164.14.00 09/16 16:51

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:52
02 Std. Sol.	0.019	16:52
03 Room Air	0.000	16:53
04 Std. Sol.	0.020	16:54
05 Room Air	0.000	16:54
06 Std. Sol.	0.019	16:55
07 Room Air	0.000	16:55

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = DR3378 Std Sol No = 201712B

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

LOW AC

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
06/18/2019

Alcohol Analyzer SN 80-004193 8164.14.00 09/16 16:57

WET CAL CHECK

Т	est		AC	Time
01	Room	Air	0.000	16:57
02	Std.	Sol.	0.293	16:58
03	Room	Air	0.000	16:59
04	Std.	Sol.	0.295	16:59
05	Room	Air	0.000	17:00
06	Std.	Sol.	0.294	17:01
07	Room	Air	0.000	17:01

 $08 \text{ Sim Temp} = 34.0^{\circ}\text{C}$

Simul Ser No = DR7089 Std Sol No = 17350

County = 08

Oper No. = 666666

Operator Signature

CHARLES EDER
HIGH AC

Remarks:

0.300 AC____ Form 106-18000

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-004193
Location = TOXL	8164.14.00 09/16
06/19/2019	09:23

Test	AC	Time	
01 Room Air	0.000	09:24	
02 *Subject Test	INT*	09:24	
03 Room Air	0.000	09:25	

*Invalid Test

Interferent Detected

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = FemaleWeight = 150

Test = DUI Cit = INTERFERENT

Dr. Lic. = ND/DIS921456

Lot No = 34917080A3

Cyl No = 7

Remarks:

Expiration Date = 02/05/2020

County = 08

Oper No. = 666666

Operator Signature ront Check +0.05% ACETONE

CHARLES EDER

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004193
Location = TOXL 8164.14.00 09/16
06/19/2019 09:25

Test		AC	Time
01 Dia	gnostic	OK	09:26
02 Room	m Air	0.000	09:27
03 *Sul	oject Test	RFI*	09:27
04 Room	n Air	0.000	09:27

*Invalid Test Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT

Sub DOB = 02/01/1992

Sub Sex = Female Weight = 150
Test = DUI Cit = RFI CHECK

Dr. Lic. = ND/DIS921456

Lot No = 34917080A3

Cyl No = 7

Expiration Date = 02/05/2020

County = 08 Oper No. = 666666

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature CHARLES EDER

Remarks:

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
06/19/2019

Alcohol Analyzer SN 80-004193 8164.14.00 09/16 09:27

DRY CAL CHECK

T	est		AC	Time
01	Room	Air	0.000	09:28
02	Std.	Gas	0.077	09:28
03	Room	Air	0.000	09:29
04	Std.	Gas	0.078	09:29
05	Room	Air	0.000	09:30
06	Std.	Gas	0.078	09:30
07	Room	Air	0.000	09:31

Lot No = 34917080A3

Cyl No = 7

Exp Date = 02/05/2020

County = 08

Oper No. = 666666

Operator Signature
CHARLES EDER

ALIBRATION CHECK

. () 1 (. .

0-080 AC Form 106-18000

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
06/19/2019

Alcohol Analyzer SN 80-004193 8164.14.00 09/16 09:31

DRY CAL CHECK

Т	est		AC	Time
01	Room	Air	0.000	09:31
02	Std.	Gas	0.078	09:32
03	Room	Air	0.000	09:32
04	Std.	Gas	0.078	09:33
05	Room	Air	0.000	09:33
06	Std.	Gas	0.078	09:34
07	Room	Air	0.000	09:34

Lot No = 34917080A3

Cyl No = 7

Exp Date = 02/05/2020

County = 08

Oper No. = 666666

ALIBRATION CAECK

Operator Signature CHARLES EDER

Remarks:

Form 106-I8000

0.080 AC

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 06/19/2019 Alcohol Analyzer SN 80-004193 8164.14.00 09/16 09:34

DRY CAL CHECK

Test		AC	Time
01 Roc	om Air	0.000	09:35
02 St	i. Gas	0.078	09:35
03 Roc	om Air	0.000	09:36
04 Sto	i. Gas	0.078	09:36
05 Roc	om Air	0.000	09:37
06 Sto	i. Gas	0.078	09:37
07 Roc	om Air	0.000	09:37

Lot No = 34917080A3

Cyl No = 7

Exp Date = 02/05/2020

County = 08

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

Form 106-I8000

0-080 AC