



**ETHANOL GAS STANDARD CYLINDER REPORT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Vanessa Lyn Richter</i>		
Location <i>GRPD</i>	Intoxilyzer® Serial Number <i>80-005360</i>	
Gas Lot Number <i>34418080A2</i>	Gas Cylinder Number <i>59</i>	Gas Expiration Date <i>02-05-2021</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).  
 A. *0.079* AC  
 B. *0.079* AC  
 C. *0.080* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>[Signature]</i>	Date <i>08-18-2019</i>
Reviewed By (Crime Laboratory Use Only) <i>[Signature]</i>	Date <i>8.26.19</i>

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005360  
Location = GFPD      8164.14.00 09/16  
08/18/2019      14:42

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:43
02 Std. Gas	0.079	14:44
03 Room Air	0.000	14:44
04 Std. Gas	0.079	14:45
05 Room Air	0.000	14:45
06 Std. Gas	0.080	14:45
07 Room Air	0.000	14:46

Lot No = 34418080A2

Cyl No = 59

Exp Date = 02/05/2021

County = 18

Oper No. = 130674



Operator Signature  
VANESSA RICHTER

Remarks: ACA For set up

Form 106-I8000