

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-005945	Instrument Location		
Reason for Install/Repair			
Install After Receiving From Crime Laboratory			
Other (Specify)			
Check When Done:			
X 1. Surge Protector Installed/Property Grounded.			
2. Telephone Line Connected to Intoxilyzer® 8000.			
∑ ⁷ 3. Breath Tube Heated.			
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).			
S. Scan/Enter Gas Cylinder Information (Level 1, Function S).			
∑ 6. Run Tests:			
A. Print Test (Level1, Function P).			
B. ACA Test (Level 1, Function C).			
🔀 C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).			
7. Repair and/or Maintenance Performed (if any):			
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.			
X 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.			
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	kilyzer® 8000 Installation and Repair C	heckout (SFN59281, Form	
Field Inspector Signature		Date	
The Ballie		7-22-20	

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

	$\gamma \Lambda$	
Reviewed/Certified By	Charles Each	Certified Date 7• 22• 2020

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005945 Location = BSPD 8164.14.00 09/16 07/22/2020 10:54 abcdefghijklmnopgrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 24119080A1 Standard Cyl #? 31 Standard Expiration? 11/05/2021 Oper No? 069823 Flow Cal. Date: 06/20/2013 Slope 637 Intercept -566407 IR Calibration Date: 05/22/2015 3um 9um ------0th Coef(*100):-21213-230931st Coef(*100):2791711392922nd Coef(*100):28981291H2O adj(mg/l*10k):332319 ********** Printer Test End **********

1- Bur

Operator Signature TIMOTHY BLETH

Remarks: " Priort Test new Install"

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-005945
Location = BSPD	8164.14.00 09/16
07/22/2020	10:54

			DRY CAL	CHECK	
Te	est			AC	Time
02 03 04 05	Room Std. Room Std. Room Std.	Gas Air Gas Air		0.000 0.078 0.000 0.078 0.000 0.079	10:55 10:55 10:56 10:56 10:57 10:57
	Room			0.000	10:57

Lot No = 24119080A1Cyl No = 31Exp Date = 11/05/2021County = 08Oper No. = 069823

T- Ble

Operator Signature TIMOTHY BLETH

Remarks: ('ACATest" New Install

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-005945 Location = BSPD 8164.14.00 09/16 07/22/2020 10:58

DRY CAL	CHECK	
Test	AC	Time
01 Room Air 02 Room Air	RFI* 0.000	11:00 11:00
*Invalid Test Inhibited - RFI		
Lot No = 24119080A1 Cyl No = 31 Exp Date = 11/05/2021 County = 08	Oper No.	= 069823

Ble

Operator Signature TIMOTHY BLETH

Remarks: "RFI Test " New install

Form 106-18000



INTOXILYZER® RECORD OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION SFN 50496 (5-2011)

Intoxilyzer® Serial	Number - 00 59 45	Location	SPPO
Gas Lot Number 2 41/9080 A1		Gas Cylinder Number 3/	Gas Expiration Date
Test Date	Chemical Test Operator Number	Operator's Name r (PRINT Last Name, First, MI)	
7-22-20	069823	Bleth, Timothy J	"Print Test, New Install"
7-22-20	069823	Heth Troiothy J.	"ACA, New Install"
7-22-20	069823	Bleth Timothy 5	"ACA, New Install" "ACA, New Install" "RFI Test, New Install"
	,		
			·

Continue on Back ->

Page 5 of 5 AGENCY RETAIN AT INDOXIDE 7822/2020 CENE