



ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Pamela B. Pehl</i>		
Location <i>Stanton - ND</i>	Intoxilyzer® Serial Number <i>80-004942</i>	
Gas Lot Number <i>05620080 A1</i>	Gas Cylinder Number <i>CYL 031</i>	Gas Expiration Date <i>04/05/2022</i>

Check When Done:

- 1 Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2 Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
A. *0.082* AC
B. *0.081* AC
C. *0.081* AC
- 3 Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4 File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5 Send the Following to the Crime Laboratory:
A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
B. ACA Test Record.

Chemical Test Operator Signature <i>P. Pehl</i>	Date <i>7/30/20</i>
Reviewed By (Crime Laboratory Use Only) <i>Charles E. Ed</i>	Date <i>7/30/20</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004942
Location = STAT 8164.14.00 09/16
07/30/2020 14:07

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:08
02 Std. Gas	0.082	14:08
03 Room Air	0.000	14:09
04 Std. Gas	0.081	14:09
05 Room Air	0.000	14:10
06 Std. Gas	0.081	14:10
07 Room Air	0.000	14:11

Lot No = 05620080A1
Cyl No = 31
Exp Date = 04/05/2022
County = 29 Oper No. = 120906



Operator Signature
PAMELA PEHL

Remarks:

Form 106-I8000