



**ETHANOL GAS STANDARD CYLINDER REPORT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Samuel A Forstre</i>		
Location <i>Wahpet City Police Department</i>	Intoxilyzer® Serial Number <i>80-004956</i>	
Gas Lot Number <i>05620080A1</i>	Gas Cylinder Number <i>CYL 001</i>	Gas Expiration Date <i>4/15/2022</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).  
 A. *0.081* AC  
 B. *0.082* AC  
 C. *0.081* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>[Signature]</i>	Date <i>8/28/2020</i>
Reviewed By (Crime Laboratory Use Only) <i>[Signature]</i>	Date <i>9.1.2020</i>


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer  
North Dakota Model 8000 SN 80-004956  
Location = WCPD 8164.14.00 09/16  
08/28/2020 19:17

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	19:18
02 Std. Gas	0.081	19:18
03 Room Air	0.000	19:19
04 Std. Gas	0.080	19:19
05 Room Air	0.000	19:20
06 Std. Gas	0.081	19:20
07 Room Air	0.000	19:21

Lot No = 05620080A1  
Cyl No = 01  
Exp Date = 04/05/2022  
County = 27 Oper No. = 131815

  
Operator Signature  
SAMUEL FORSTIE

Remarks: Good Test

Form 106-I8000