



ETHANOL GAS STANDARD CYLINDER REPORT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Kathryn Allen</i>		
Location <i>Cass County Sheriff's Office - Jail</i>	Intoxilyzer® Serial Number <i>80-006496</i>	
Gas Lot Number <i>0562080A1</i>	Gas Cylinder Number <i>27</i>	Gas Expiration Date <i>04/05/2022</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 A. 0. 082 AC
 B. 0. 082 AC
 C. 0. 082 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>[Signature]</i> # <i>3071</i>	Date <i>04-14-2021</i>
Reviewed By (Crime Laboratory Use Only) <i>[Signature]</i>	Date <i>4.26.21</i>


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006496
Location = CASS 8164.13.00 06/09
04/14/2021 18:35

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	18:36
02 Std. Gas	0.082	18:36
03 Room Air	0.000	18:37
04 Std. Gas	0.082	18:37
05 Room Air	0.000	18:38
06 Std. Gas	0.082	18:38
07 Room Air	0.000	18:38

Lot No = 05620080A1
Cyl No = 27
Exp Date = 04/05/2022
County = 09 Oper No. = 090801


Operator Signature
KATHRYN ALLEN

Remarks:

Form 106-I8000