



ETHANOL GAS STANDARD CYLINDER REPORT
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
 SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Bleth, Timothy J</i>		
Location <i>BSPD</i>	Intoxilyzer® Serial Number <i>80-005945</i>	
Gas Lot Number <i>02621080 A1</i>	Gas Cylinder Number <i>20</i>	Gas Expiration Date <i>3-5-23</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C). *RECEIVED ACA TEST. CEE*
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
CEE 5-6-21
 A. *0.080* AC
 B. *0.079* AC
 C. *0.079* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>T. Bleth</i>	Date <i>5-6-21</i>
Reviewed By (Crime Laboratory Use Only) <i>Charles E. Ede</i>	Date <i>5-6-21</i>


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005945
Location = BSPD 8164.14.00 09/16
05/06/2021 13:07

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:08
02 Std. Gas	0.080	13:08
03 Room Air	0.000	13:09
04 Std. Gas	0.079	13:09
05 Room Air	0.000	13:10
06 Std. Gas	0.079	13:10
07 Room Air	0.000	13:11

Lot No = 02621080A1
Cyl No = 20
Exp Date = 03/05/2023
County = 08 Oper No. = 069823


Operator Signature
TIMOTHY BLETH

Remarks: *new install & cylinder change "ACA"*

Form 106-I8000