

ETHANOL GAS STANDARD CYLINDER REPORT NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
SFN 59282 (1987)

Chemical Test Operator Name (Print)			
Location			
Fargo Police PD 3 1230 2	1.5th St. 5. Fyo 80-00	D649 <i>5</i>	
02621080A1	Gas Cylinder Number	Gas Expiration Date (03/05/2023	3

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C). Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0.081 AC
 - B. 0.082 AC
 - C. 0.082AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature	
7711	Date
Fette Jane	05/30/2021
Reviewed By (Crime Laboratory Use Only)	
land I	Date 7/22/
Julie to the	1122121
	1 /01

FORM 105-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = FAPD
05/30/2021

Alcohol Analyzer SN 80-006495 8164.14.00 09/16 22:45

DRY CAL CHECK

Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas	0.000 0.081 0.000 0.082 0.000 0.082 0.000	22:46 22:46 22:47 22:47 22:48 22:48 22:49
07 Room Air	0.000	

Lot No = 02621080A1

Cyl No = 13

Exp Date = 03/05/2023

County = 09

Oper No. = 130056

Operator Signature PATRICK THOMAS

Remarks: ACA Calibration Test for Installation.

Form 106-I8000