



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

*DUNN aew 20Jan23*

Serial Number <b>80-005955</b>	Instrument Location <b>DCSO</b>
Reason for Install/Repair	
<input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input type="checkbox"/> Other (Specify) _____	

*Dunn County Sheriff's Office - ~~Dunn~~ however there is already a DCSO elsewhere*

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): None
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature 	Date <b>06-21-2021</b>
-------------------------------	---------------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By <i>Janelle Portscheller</i>	Certified Date <b>20 Apr 23</b>
--	------------------------------------

Intoxilyzer Test Record and Checklist  
 NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer           Alcohol Analyzer  
 North Dakota Model 8000       SN 80-005955  
 Location = DUNN               8164.14.00 09/16  
 06/21/2021                               11:29

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-|=|  
 ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-|=|  
 ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#%&^\*()\_+?

Current Instrument Setup

Data Entry Mode:               Enabled  
 Start Test Sequence:         DABACABA  
 Display Prelim Rslt?         Yes  
 Display Third Digit?         Yes  
 Inhib Printer(Y/N)?         No  
 Display Volume?               No  
 Disable On Memfull?         Yes  
 # of Print Copies?           1  
 Select Std (D/W/I)?         Dry  
 Standard Value?               0.080  
 Standard Lot #?               05620080A1  
 Standard Cyl #?               35  
 Standard Expiration?         04/05/2022  
 Oper No?                       132823

Flow Cal. Date:               06/21/2013  
 Slope                         757  
 Intercept                     -837385

IR Calibration Date:         06/21/2013  
                                   3um                       9um

-----  
 0th Coef(\*100):         -19897               -22878  
 1st Coef(\*100):         268381               131979  
 2nd Coef(\*100):         5157                   1473  
 H2O adj(mg/l\*10k):    717                     600

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
 EASTON BROST

Remarks:

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-005955  
Location = DUNN                8164.14.00 09/16  
06/21/2021                        11:35

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:36
02 Std. Gas	0.083	11:36
03 Room Air	0.000	11:36
04 Std. Gas	0.083	11:37
05 Room Air	0.000	11:37
06 Std. Gas	0.082	11:38
07 Room Air	0.000	11:38

Lot No = 05620080A1  
Cyl No = 35  
Exp Date = 04/05/2022  
County = 13                            Oper No. = 132823

  
Operator Signature  
EASTON BROST

Remarks:

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005955  
Location = DUNN      8164.14.00 09/16  
06/21/2021      11:30

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:30
02 Std. Gas	0.082	11:31
03 Room Air	0.000	11:31
04 Std. Gas	0.082	11:32
05 Room Air	0.000	11:32
06 Std. Gas	0.082	11:33
07 Room Air	RFI*	11:33
08 Room Air	0.000	11:33

\*Invalid Test  
Inhibited - RFI

Lot No = 05620080A1  
Cyl No = 35  
Exp Date = 04/05/2022  
County = 13      Oper No. = 132823



Operator Signature  
EASTON BROST

Remarks: