

ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59282 (08/2018)

Chemical Test Operator Name (Print) Jared Lemieux					
		Intoxilyzer® Serial Number 80-005358			
Gas Lot Number 9	Gas Cylinder Number 02621080A1		Gas Expiration Date 3/5/2023		

Check When Done:

- X 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).

Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).

- A. 0.079 AC
- B. 0.080 AC
- C. 0.079 AC
- 🗵 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- X 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.

	5. Send the Following to the Crime Laboratory:	RECENED	FORMS & TEST RECORDS
	A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).		

A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G). B. ACA Test Record.

Chemical Test Operator Signature 4834	Date 9/30/2021
Reviewed By (Crime Laboratory Use Only)	Date 9/30/2021

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000Alcohol Analyzer
SN 80-005358
8164.14.00 09/16
15:50

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.079 0.000 0.080 0.000 0.079 0.000	15:51 15:52 15:52 15:53 15:53 15:53 15:53

Lot No = 02621080A1 Cyl No = 9 Exp Date = 03/05/2023 County = 08 · Oper No. = 130315

1834 Operator Signature

JARED LEMIEUX

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 Alcohol Analyzer CMI, Inc. Intoxilyzer North Dakota Model 8000 SN 80-005358 8164.14.00 09/16 Location = BMDC15:48 09/27/2021 *********** Printer Test **************** abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPORSTUVWXYZ!@#\$%^&*()_+? abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+? Current Instrument Setup Enabled Data Entry Mode: DABACABA Start Test Sequence: Display Prelim Rslt? Yes Yes Display Third Digit? No Inhib Printer(Y/N)? No Display Volume? Yes Disable On Memfull? 1 # of Print Copies? Select Std (D/W/I)? Dry 0.080 Standard Value? 02621080A1 Standard Lot #? 9 Standard Cyl #? 03/05/2023 Standard Expiration? 130315 Oper No? 06/04/2019 Flow Cal. Date: 570 Slope 4079 Intercept 06/04/2019 IR Calibration Date: 9um 3um -------24652 -19471 0th Coef(*100): 143582 1st Coef(*100): 273237 1293 2965 2nd Coef(*100): 374 H2O adj(mg/l*10k): 273

********** Printer Test End **********

Operator Signature

JARED LEMIEUX

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000Alcohol Analyzer
SN 80-005358Location = BMDC8164.14.00 09/16
15:57

DRY CA	L CHECK	
Test	AC	Time
01 Room Air 02 Room Air	RFI* 0.000	15:58 15:58
*Invalid Test Inhibited - RFI		
Lot No = 02621080A1 Cyl No = 9		

Cyl No = 9 Exp Date = 03/05/2023 County = 08 Oper No. = 130315

1834 Operator Signature

JARED LEMIEUX

Remarks:

Form 106-18000