



ETHANOL GAS STANDARD CYLINDER REPORT
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
 SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Dylan Thomas 4402</i>		
Location <i>ASHL</i>	Intoxilyzer® Serial Number <i>80-006490</i>	
Gas Lot Number <i>0262/080A1</i>	Gas Cylinder Number <i>36</i>	Gas Expiration Date <i>03/05/2023</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).

A. 0.____ AC	<i>0.079 AC</i>	<i>CEE</i>
B. 0.____ AC	<i>0.079 AC</i>	<i>11.2.21</i>
C. 0.____ AC	<i>0.079 AC</i>	<i>RECEIVED ACA TEST (DRY CAL CHECK)</i>
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Dylan Thomas 4402</i>	Date <i>10/22/2021</i>
Reviewed By (Crime Laboratory Use Only) <i>Charles E. Ehl</i>	Date <i>11/2/21</i>


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006490
Location = ASHL 8164.14.00 09/16
10/28/2021 12:51

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:52
02 Std. Gas	0.079	12:52
03 Room Air	0.000	12:53
04 Std. Gas	0.079	12:53
05 Room Air	0.000	12:54
06 Std. Gas	0.079	12:54
07 Room Air	0.000	12:55

Lot No = 02621080A1
Cyl No = 36
Exp Date = 03/05/2023
County = 26 Oper No. = 132584


Operator Signature
DYLAN THOMES

Remarks:

Form 106-I8000

THIS IS A TRUE AND ACCURATE COPY.

