

ETHANOL GAS STANDARD CYLINDER REPORT NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59282 (08/2018)

| Chemical Test Operator Name (Print) Jared Lemieux | | | | | |
|--|---------------------|----------------------------|---------------------|--|--|
| Location | | Intoxilyzer® Serial Number | | | |
| BMDC | | 80-004945 | | | |
| Gas Lot Number | Gas Cylinder Number | | Gas Expiration Date | | |
| 02621080A1 | 7 | | 03/05/2023 | | |

Check When Done:

- X 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- X 2. Perform an ACA Test (Level 1, Function C).

Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).

- A. 0.06L AC
- B. 0.28 L AC
- C. 0.08 0 AC
- X 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- X 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- X 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

| Chemical Test Operator Signature 4834 | Date 10/28/2021 |
|---|--------------------|
| Reviewed By (Crime Laboratory Use Only) | Date 11.2.2021 |

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-004945Location = BMDC8164.14.00 09/1610/28/202109:59

| DRY CAL CHECK | | | | |
|---|---|---|--|--|
| Test | AC | Time | | |
| 01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air | 0.000 0.081 0.000 0.081 0.000 0.080 0.000 | 10:00 10:00 10:01 10:01 10:01 10:02 10:02 | | |

Lot No = 02621080A1 Cyl No = 7 Exp Date = 03/05/2023 County = 08 Oper No. = 130315

1834 Operator Signature

JARED LEMIEUX

Remarks:

Form 106-18000