

## ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59282 (08/2018)

Chemical Test Operator Name (Print) DELBERT HOFFMANN					
Location SHER		Intoxilyzer® Serial Number 20 - 005950			
Gas Lot Number 260 2 1 0 8 0 A 1	Gas Cylinder Number	34	Gas Expiration Date		

## Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C). Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
  - A. 0.082 AC B. 0.082 AC C. 0.081 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature		Date 11-04- 2021
Reviewed By (Crime Laberlatory Use Only)	Clearly Est	Date 11. 4. 2021

FORM 105-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = SHER
11/04/2021

Alcohol Analyzer SN 80-005950 8164.14.00 09/16 08:44

## DRY CAL CHECK

Test	AC	Time	
01 Room Air	0.000	08:45	
02 Std. Gas	0.082	08:45	
03 Room Air	0.000	08:46	
04 Std. Gas	0.082	08:46	
05 Room Air	0.000	08:46	
06 Std. Gas	0.081	08:47	
07 Room Air	0.000	08:47	

Lot No = 26021080A1

Cyl No = 34

Exp Date = 10/05/2023

County = 42

Oper No. = 048909

Operato Signature DELBERT HOFFMANN

Remarks:

ACA TEST

Form 106-I8000