



**ETHANOL GAS STANDARD CYLINDER REPORT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59282 (08/2018)

|   |   |                                       |
|---|---|---------------------------------------|
| Chemical Test Operator Name (Print) <i>Gabriel Johnston</i> |   |                                       |
| Location <i>Candler County</i>                              | Intoxilyzer® Serial Number <i>80-005951</i> |                                       |
| Gas Lot Number <i>26021080A1</i>                            | Gas Cylinder Number <i>29</i>               | Gas Expiration Date <i>10-05-2023</i> |

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).  
 A. 0. *062* AC  
 B. 0. *083* AC  
 C. 0. *082* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

|   |                        |
|---|------------------------|
| Chemical Test Operator Signature <i>[Signature]</i> # <i>4228</i> | Date <i>11-09-2021</i> |
| Reviewed By (Crime Laboratory Use Only) <i>[Signature]</i>        | Date <i>11/16/2021</i> |


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005951  
Location = LANG      8164.14.00 09/16  
11/09/2021      18:55

DRY CAL CHECK

| Test        | AC    | Time  |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 18:56 |
| 02 Std. Gas | 0.082 | 18:56 |
| 03 Room Air | 0.000 | 18:57 |
| 04 Std. Gas | 0.083 | 18:57 |
| 05 Room Air | 0.000 | 18:58 |
| 06 Std. Gas | 0.082 | 18:58 |
| 07 Room Air | 0.000 | 18:59 |

Lot No = 26021080A1  
Cyl No = 29  
Exp Date = 10/05/2023  
County = 10      Oper No. = 021009

  
\_\_\_\_\_  
Operator Signature  
GABRIEL JOHNSTON

Remarks:

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005951  
Location = LANG      8164.14.00 09/16  
11/09/2021      18:51

DRY CAL CHECK

| Test        | AC    | Time  |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 18:51 |
| 02 Std. Gas | 0.000 | 18:52 |
| 03 Room Air | 0.000 | 18:52 |
| 04 Std. Gas | 0.083 | 18:53 |
| 05 Room Air | 0.000 | 18:53 |
| 06 Std. Gas | 0.082 | 18:53 |
| 07 Room Air | 0.000 | 18:54 |

Lot No = 26021080A1  
Cyl No = 29  
Exp Date = 10/05/2023  
County = 10      Oper No. = 021009

  
\_\_\_\_\_  
Operator Signature  
GABRIEL JOHNSTON

Remarks: *started before putting tank on. (G) 4228*

Form 106-I8000