

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-005358		Instrument Location BMDC		
Reason for Install/Repair				
Install After Receiving From Crime Laboratory [Install Afte	er Location Change		

X	Other	(Specify)	December	Routine	Checks
---	-------	-----------	----------	---------	--------

Check When Done:

- X 1. Surge Protector Installed/Property Grounded.
- X 2. Telephone Line Connected to Intoxilyzer® 8000.
- X 3. Breath Tube Heated.
- X 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- X 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).

X 6. Run Tests:

- X A. Print Test (Level1, Function P).
- X B. ACA Test (Level 1, Function C).
- C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- X 7. Repair and/or Maintenance Performed (if any):
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- X 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature 4834	Date 12/5/2021

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By	11/1/1	Certified	Pate	1	
	Maderta	12/	23/	2021	
		7		a manufacture of the second states and the second states and	

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer SN 80-005358 North Dakota Model 8000 8164.14.00 09/16 Location = BMDC00:59 12/05/2021 abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+? abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+? Current Instrument Setup Enabled Data Entry Mode: DABACABA Start Test Sequence: Yes Display Prelim Rslt? Display Third Digit? Yes No Inhib Printer(Y/N)? No Display Volume? Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry 0.080 Standard Value? 02621080A1 Standard Lot #? Standard Cyl #? 9 03/05/2023 Standard Expiration? 130315 Oper No? 06/04/2019 Flow Cal. Date: 570 Slope 4079 Intercept 06/04/2019 IR Calibration Date: 9um 3um -24652 -19471 0th Coef(*100): 2nd Coef(*100): H2O add (- (*100): 143582 273237 1293 2965 H2O adj(mg/l*10k): 273 374

********* Printer Test End ***********

Operator Signature JARED LEMIEUX

Remarks:

December Checks (2021)

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-005358Location = BMDC8164.14.00 09/1612/05/202101:00

			DRY CAL CHECK]
$T\epsilon$	est		AC	Time
03 04 05 06	Room Std. Room Std. Room Std. Room	Gas Air Gas Air Gas	0.000 0.082 0.000 0.081 0.000 0.082 0.000	01:01 01:01 01:02 01:02 01:02 01:03 01:03

Lot No = 02621080A1 Cyl No = 9 Exp Date = 03/05/2023 County = 08 Oper No. = 130315

934 Operator Stonature

JARED LEMIEUX

Remarks:

rks: Dre. Routine chuk (2021)

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000Alcohol Analyzer
SN 80-005358Location = BMDC8164.14.00 09/1612/05/202101:07

DRY CA	L CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Room Air	0.000 0.082 RFI* RFI*	01:07 01:07 01:08 01:08
*Invalid Test Inhibited - RFI		
Lot No = 02621080A1 Cyl No = 9		

 Cyr No = 9

 Exp Date = 03/05/2023

 County = 08

 Oper No. = 130315

4834 Operator Ire Signat JARED LEMIEUX

Remarks:

Dre. Routine charle (2021)

Form 106-18000