

ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59282 (08/2018)

Chemical Test Operator Name (Mark New Holy	Print)		
Location Three Affiliated Tribes- New Town		Intoxilyzer® Serial Number 80-005940	
Gas Lot Number 26021080A1	Gas Cylinder Number 15	Gas Expiration Date 10/05/2023	

Check When Done:

- X 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C). Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0.0 78 AC
 - B. 0.2 3 \$ AC
 - C. 0.0 77 AC
- 🗵 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- S. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature	
Reviewed By (Crime Laboratory Use Only)	Date 7/27/2022
Jarelle Portschuler	Date 2 34 / 3023

* Note Form 105-G was filled out by Mark New Holy because Shannon Dursen forgot to fill out this form and is on leave from his agency.

62/24/2023 ggs

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = NEWT
07/27/2022

Alcohol Analyzer SN 80-005940 8164.14.00 09/16 11:25

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:27
02 Std. Gas	0.078	11:27
03 Room Air	0.000	11:27
04 Std. Gas	0.078	11:28
05 Room Air	0.000	11:28
06 Std. Gas	0.078	11:29
07 Room Air	0.000	11:29

Lot No = 26021080A1

Cyl No = 15

Exp Date = 10/05/2023

County = 60

Oper No. = 020632

Operator Signature SHANNON DUISEN

Remarks:

Form 106-I8000