



ETHANOL BREATH STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (5-2011)

Chemical Test Operator Name (Print) <i>Joshua T Trevena</i>		
Location <i>Wataford City Police Dept</i>	Intoxilyzer® Serial Number <i>80-005958</i>	
Gas Lot Number <i>19615080 AS</i>	Gas Cylinder Number <i>14</i>	Gas Expiration Date <i>09.05.2017</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. *0.078* AC
 - B. *0.078* AC
 - C. *0.078* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Breath Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Joshua T Trevena 1106</i>	Date <i>07-03-2016</i> ✓
Reviewed By (Crime Laboratory Use Only) <i>Deb Shaver</i>	Date <i>05 July 2016</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005958
Location = WCPD 8164.13.00 06/09
07/03/2016 16:03

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:04
02 Std. Gas	0.078	16:04
03 Room Air	0.000	16:05
04 Std. Gas	0.078	16:05
05 Room Air	0.000	16:06
06 Std. Gas	0.078	16:06
07 Room Air	0.000	16:06

Lot No = 19615080A5

Cyl No = 14

Exp Date = 09/05/2017

County = 27

Oper No. = 130612


Operator Signature
JOSHUA TREVENA

Remarks: Calibration check after tank exchange.
Load Test
New Gas Standard Cylinder
Form 106-I8000