



ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (10-2013)

Chemical Test Operator Name (Print) David Feland		
Location Steele	Intoxilyzer® Serial Number 80-005937	
Gas Lot Number 19615080 AS	Gas Cylinder Number 27	Gas Expiration Date 09/05/2017

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0.079 AC
 - B. 0.081 AC
 - C. 0.080 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature David Feland	Date 11 Apr 16
Reviewed By (Crime Laboratory Use Only) Deb Shanaver	Date 11 Apr 2016


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005937
Location = STEE 8164.13.00 06/09
04/11/2016 13:30

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:31
02 Std. Gas	0.079	13:32
03 Room Air	0.000	13:32
04 Std. Gas	0.081	13:32
05 Room Air	0.000	13:33
06 Std. Gas	0.080	13:33
07 Room Air	0.000	13:34

Lot No = 19615080A5
Cyl No = 27
Exp Date = 09/05/2017
County = 22 Oper No. = 119404


Operator Signature
DAVID FELAND

Remarks: