



# ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION  
Toxicology Section/Breath Alcohol Program  
SFN 59282 (10-2013)

Chemical Test Operator Name (Print) David Feland		
Location Kidder J.O.	Intoxilyzer® Serial Number 80-005951	
Gas Lot Number 33913080A4	Gas Cylinder Number 48	Gas Expiration Date 01/01/16

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).  
A. 0.081 AC  
B. 0.081 AC  
C. 0.080 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature David Feland	Date 15 Oct 14
Reviewed By (Crime Laboratory Use Only) DJP Deb Kashner	Date 17 Oct 14

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005951  
Location = STEE      8164.13.00 06/09  
10/15/2014      13:23

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:25
02 Std. Gas	0.081	13:25
03 Room Air	0.000	13:26
04 Std. Gas	0.081	13:26
05 Room Air	0.000	13:27
06 Std. Gas	0.080	13:27
07 Room Air	0.000	13:28

Lot No = 33913080A4  
Cyl No = 48  
Exp Date = 01/01/2016  
County = 22      Oper No. = 119404



Operator Signature  
DAVID FELAND

Remarks:

Form 106-I8000