



# ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL



CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <b>JON WOLD.</b>		
Location <b>SARGENT county Sheriff's office.</b>	Intoxilyzer® Serial Number <b>80-000690</b>	
Gas Lot Number <b>34921080 A1</b>	Gas Cylinder Number <b>57</b>	Gas Expiration Date <b>02/05/</b>

Check When Done.

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).  
A. 0.080 AC  
B. 0.080 AC  
C. 0.080 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G)
  - B. ACA Test Record.

Chemical Test Operator Signature 	Date <b>4/26/23</b>
Reviewed By (Crime Laboratory Use Only) 	Date <b>26 Apr 2023</b>

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006690  
Location = SACO      8164.14.00 09/16  
04/26/2023      16:11

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:11
02 Std. Gas	0.080	16:12
03 Room Air	0.000	16:12
04 Std. Gas	0.080	16:13
05 Room Air	0.000	16:13
06 Std. Gas	0.080	16:14
07 Room Air	0.000	16:14

Lot No = 34921080A1  
Cyl No = 57  
Exp Date = 02/05/2024  
County = 41      Oper No. = 132607



Operator Signature  
JONATHAN WOLD

Remarks:

Form 106-I8000