





**ETHANOL GAS STANDARD CYLINDER REPORT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <b>HUARD, CHRISTOPHER (132570)</b>		
Location <b>RICHLAND CO. JAIL; WAHPETON</b>	Intoxilyzer® Serial Number <b>80-006508</b>	
Gas Lot Number <b>349Z1080A1</b>	Gas Cylinder Number <b>035</b>	Gas Expiration Date <b>02/05/2024</b>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).  
 A. **0.079** AC  
 B. **0.079** AC  
 C. **0.078** AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature 	Date <b>06/06/2023</b>
Reviewed By (Crime Laboratory Use Only) 	Date <b>12 June 23</b>


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006508  
Location = WAHP      8164.14.00 09/16  
06/06/2023      16:56

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:56
02 Std. Gas	0.079	16:57
03 Room Air	0.000	16:57
04 Std. Gas	0.079	16:58
05 Room Air	0.000	16:58
06 Std. Gas	0.078	16:58
07 Room Air	0.000	16:59

Lot No = 34921080A1  
Cyl No = 35  
Exp Date = 02/05/2024  
County = 39      Oper No. = 132570

  
Operator Signature  
CHRISTOPHER HUARD

Remarks: *ACA TEST*

Form 106-I8000