



**ETHANOL GAS STANDARD CYLINDER REPORT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Tracy L. Meidinger</i>		
Location <i>UNDPD</i>	Intoxilyzer® Serial Number <i>80-005955</i>	
Gas Lot Number <i>34921080A1</i>	Gas Cylinder Number <i>30</i>	Gas Expiration Date <i>02-05-2024</i>

Check When Done:

- TM*  1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- TM*  2. Perform an ACA Test (Level 1, Function C).  
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).  
 A. *0.079* AC  
 B. *0.076* AC  
 C. *0.079* AC
- TM*  3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- TM*  4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- TM*  5. Send the Following to the Crime Laboratory:  
 A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).  
 B. ACA Test Record.

Chemical Test Operator Signature <i>Tracy L. Meidinger</i>	Date <i>6-30-23</i>
Reviewed By (Crime Laboratory Use Only) <i>Aime Nahl</i>	Date <i>03Jul2023</i>


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005955  
Location = UNDP      8164.14.00 09/16  
06/30/2023      12:46

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:47
02 Std. Gas	0.079	12:48
03 Room Air	0.000	12:48
04 Std. Gas	0.078	12:48
05 Room Air	0.000	12:49
06 Std. Gas	0.079	12:49
07 Room Air	0.000	12:50

Lot No = 34921080A1  
Cyl No = 30  
Exp Date = 02/05/2024  
County = 18      Oper No. = 119409

  
Operator Signature  
TRACY MEIDINGER

Remarks:

Form 106-I8000