



**ETHANOL GAS STANDARD CYLINDER REPORT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Jesus Armado Barrera Iracheta</i>		
Location <i>McIntosh County Correctional Center</i>	Intoxilyzer® Serial Number <i>80 - 006676</i>	
Gas Lot Number <i>01923080A3</i>	Gas Cylinder Number <i>22</i>	Gas Expiration Date <i>02-05-2025</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).  
 A. 0.080 AC  
 B. 0.080 AC  
 C. 0.082 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:  
 A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).  
 B. ACA Test Record.

Chemical Test Operator Signature <i>Jesus Armado Barrera Iracheta</i>	Date <i>04/26/23</i>
Reviewed By (Crime Laboratory Use Only) <i>Amir</i>	Date <i>27Apr2023</i>


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006676  
Location = MOSO      8164.14.00 09/16  
04/26/2023      09:50

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:51
02 Std. Gas	0.080	09:51
03 Room Air	0.000	09:52
04 Std. Gas	0.080	09:52
05 Room Air	0.000	09:53
06 Std. Gas	0.082	09:53
07 Room Air	0.000	09:53

Lot No = 01923080A3  
Cyl No = 22  
Exp Date = 02/05/2025  
County = 31      Oper No. = 133127



Operator Signature  
JESUS BARRERA

Remarks:

Form 106-I8000