

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number Q.o. Instrument Location	
0-006676 Mantil Charles	in Center
Reason for Install/Repair	ona Centr
Install After Receiving From Crime Laboratory Install After Location Change	
Dother (Specify) Annual Inspection at Chine Laboratory_	
Check When Done:	
1. Surge Protector Installed/Property Grounded.	
1√2. Telephone Line Connected to Intoxilyzer® 8000.	
ろ. Breath Tube Heated.	
A. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).	
X 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).	
(y) 6. Run Tests:	
A. Print Test (Level1, Function P).	
B. ACA Test (Level 1, Function C).	
C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio Du	uring Test)
7. Repair and/or Maintenance Performed (if any):	
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the I	ntoxilvzer® for Use
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Age	
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair C 104-G), Print Test, ACA Test, and RFI Test.	
Field Inspector Signature	Date
tsist	4/26/23
Crime Laboratory Use Only	
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to deter the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crin County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.	rmine alcohol concentration from ne Laboratory Division, in the
Reviewed/Certified By	Certified Date

Jesus	Wasny	able to	complete	RFI	test. Shayla
Co	mpleted	test. AE	N 27Apra	123	

FORM 104-G

27Apra023

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006676 Location = MOSO 8164.14.00 09/16 09:45 04/26/2023 abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +? abcdefghijklmnopqrstuvwxyz1234567890-= ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+? Current Instrument Setup Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Display Third Digit? Inhib Printer(Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? Standard Lot #? Standard Cyl #? 0.080 01923080A3 22 02/05/2025 Standard Expiration? Oper No? 133127 08/22/2016 Flow Cal. Date: 710 Slope Intercept -713681 IR Calibration Date: 08/22/2016 3um 9um _____ 0th Coef(*100):-166671st Coef(*100):2591902nd Coef(*100):3393 -22327 132291 1374 546 H2O adj(mg/l*10k): 611

********** Printer Test End **********

Operator Signature JESUS BARRERA

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. IntoxilyzerAlcohol AnalyzerNorth Dakota Model 8000SN 80-006676Location = MOSO8164.14.00 09/1604/26/202309:50

	DRY CAL CHECK	
Test	AC	Time
01 Room Air 02 Std. Gas 03 Room Air 04 Std. Gas 05 Room Air 06 Std. Gas 07 Room Air	0.000 0.080 0.000 0.080 0.000 0.082 0.000	09:51 09:52 09:52 09:53 09:53 09:53

Lot No = 01923080A3 Cyl No = 22 Exp Date = 02/05/2025 County = 31 Oper No. = 133127

Operator Signature JESUS BARRERA

Remarks:

Form 106-18000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501 CMI, Inc. Intoxilyzer Alcohol Analyzer North Dakota Model 8000 SN 80-006676 Location = MOSO 8164.14.00 09/16 12:42 04/27/2023 Time AC Test 12:44 0.000 01 Room Air 12:45 02 *Subject Test RFI* 12:45 03 Room Air 0.000 *Invalid Test Inhibited - RFI Sub Name = TEST, TEST TEST Sub DOB = 08/19/2001Weight = 112Sub Sex = Male Cit = 11111 Test = DUI Dr. Lic. = ND/TES011222Lot No = 01923080A3Cyl No = 22Expiration Date = 02/05/2025Oper No. = 132957County = 31

1477 Signature Op rator SHAYLA RODGERS

Remarks:

Form 106-18000