



ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

Chemical Test Operator Name (Print) FRANK EBERLE		
Location Bowman County Sheriff's Office	Intoxilyzer® Serial Number 80-003071	
Gas Lot Number 14323080A1	Gas Cylinder Number 14	Gas Expiration Date 6-5-2025

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
A. 0.079 AC
B. 0.079 AC
C. 0.078 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
B. ACA Test Record.

Chemical Test Operator Signature <i>Frank Eberle</i>	Date 9-21-2023
Reviewed By (Crime Laboratory Use Only) <i>Aime Nahl</i>	Date 22Sep2023


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003071
Location = BOWM 8164.14.00 09/16
09/21/2023 14:01

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:02
02 Std. Gas	0.079	14:02
03 Room Air	0.000	14:03
04 Std. Gas	0.079	14:03
05 Room Air	0.000	14:04
06 Std. Gas	0.078	14:04
07 Room Air	0.000	14:04

Lot No = 14323080A1
Cyl No = 14
Exp Date = 06/05/2025
County = 06 Oper No. = 130015



Operator Signature
FRANK EBERLE

Remarks:

Form 106-I8000