

ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59282 (08/2018)

Chemical Test Operator Name (Print)	i Stuart			
Location 1211 12th St SE, Suite C. Watford City, ND 50054 Intoxilyzer® Serial Number 200-005362				
Gas Lot Number 4323000 A4	Gas Cylinder Number \G	Gas Expiration Date U 5 2025		

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C). Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0.0 0 4 AC B. 0.0 0 7 AC C. 0.0 0 4 AC
- 🗓 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record

	all
Chemical Test Operator Signature	Date 132324
Reviewed By (Crime Laboratory Use Only)	Date 1/17/2024

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = MCSO 01/03/2024

Alcohol Analyzer SN 80-005362 8164.14.00 09/16 07:45

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	07:46
02 Std. Gas	0.084	07:46
03 Room Air	0.000	07:47
04 Std. Gas	0.083	07:47
05 Room Air	0.000	07:47
06 Std. Gas	0.084	07:48
07 Room Air	0.000	07:48

Lot No = 14323080A4

Cyl No = 15

Exp Date = 06/05/2025

County = 27

Oper No. = 132007

Operator Signature KARI STUART

Remarks:

STATE OF NORTH DAKOTA $\}$ ss COUNTY OF MCKENZIE

I certify that the foregoing is a true and correct copy of the original document on file in this office. Dated at Watford City. North Dakota

Form 106-I8000