



ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Joshua Gudwanger</i>		
Location <i>Ward city Jail</i>	Intoxilyzer® Serial Number <i>98-006493</i>	
Gas Lot Number <i>44</i>	Gas Cylinder Number <i>H3 23080A1</i>	Gas Expiration Date <i>06-05-2025</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0.082 AC
 - B. 0.083 AC
 - C. 0.082 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>[Signature]</i> 252	Date <i>01/06/2024</i>
Reviewed By (Crime Laboratory Use Only) <i>[Signature]</i>	Date <i>08Jan2024</i>


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006493
Location = WARD 8164.14.00 09/16
01/06/2024 18:43

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	18:44
02 Std. Gas	0.083	18:44
03 Room Air	0.000	18:45
04 Std. Gas	0.083	18:45
05 Room Air	0.000	18:46
06 Std. Gas	0.082	18:46
07 Room Air	0.000	18:47

Lot No = 14323080A4
Cyl No = 44
Exp Date = 06/05/2025
County = 51 Oper No. = 131334



Operator Signature
JOSHUA GUDVANGEN

Remarks:

Form 106-I8000