



ETHANOL GAS STANDARD CYLINDER REPORT
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
 SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Jacob Nikko</i>		
Location <i>Grafton PD</i>	Intoxilyzer® Serial Number <i>80-003067</i>	
Gas Lot Number <i>01923080A3</i>	Gas Cylinder Number <i>16</i>	Gas Expiration Date <i>02/05/2025</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 A. *0.080* AC
 B. *0.080* AC
 C. *0.080* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Jacob Nikko</i>	Date <i>1/26/2025</i>
Reviewed By (Crime Laboratory Use Only) <i>Aime Nahl</i>	Date 29Jan2024


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003067
Location = GFTN 8164.14.00 09/16
01/26/2024 12:05

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:05
02 Std. Gas	0.080	12:06
03 Room Air	0.000	12:06
04 Std. Gas	0.080	12:07
05 Room Air	0.000	12:07
06 Std. Gas	0.080	12:08
07 Room Air	0.000	12:08

Lot No = 01923080A3
Cyl No = 16
Exp Date = 02/05/2025
County = 50 Oper No. = 132597



Operator Signature
JACOB NIKKO

Remarks:

Form 106-I8000

The undersigned, having custody of the original public record, certifies that the information hereon is a true and correct copy of the original document maintained as a part of the activities of this agency.

Signature/Title  / Police Officer

Agency/Date Grafton PD / 1/25/24