



**ETHANOL GAS STANDARD CYLINDER REPORT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59282 (08/2018)

Chemical Test Operator Name (Print) Jeffrey Bauske		
Location Rolla	Intoxilyzer® Serial Number 80-007097	
Gas Lot Number 14323080A4	Gas Cylinder Number 10	Gas Expiration Date 06/05/2025

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
  - A. 0.083 AC
  - B. 0.083 AC
  - C. 0.082 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>Jeffrey Bauske 287</i>	Date 2/4/2024
Reviewed By (Crime Laboratory Use Only) <i>Aimee Nahl</i>	Date 05Feb2024


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-007097  
Location = ROLL      8164.16.00 09/18  
02/04/2024      13:26

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:27
02 Std. Gas	0.083	13:27
03 Room Air	0.000	13:27
04 Std. Gas	0.083	13:28
05 Room Air	0.000	13:28
06 Std. Gas	0.082	13:29
07 Room Air	0.000	13:29

Lot No = 14323080A4  
Cyl No = 10  
Exp Date = 06/05/2025  
County = 40      Oper No. = 120501

  
Operator Signature  
JEFFREY BAUSKE

Remarks: *ACA TEST / NEW CYLINDER*

Form 106-I8000