



ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Emily Henley</i>		
Location <i>Minot AFB</i>	Intoxilyzer® Serial Number <i>80-006501</i>	
Gas Lot Number <i>14323080A4</i>	Gas Cylinder Number <i>5</i>	Gas Expiration Date <i>06-05-2025</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. *0.080* AC
 - B. *0.079* AC
 - C. *0.080* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>[Signature]</i>	Date <i>08 Feb 2024</i>
Reviewed By (Crime Laboratory Use Only) <i>[Signature]</i>	Date <i>08Feb2024</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006501
Location = MNAF 8164.14.00 09/16
02/08/2024 11:03

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:04
02 Std. Gas	0.080	11:04
03 Room Air	0.000	11:04
04 Std. Gas	0.079	11:05
05 Room Air	0.000	11:05
06 Std. Gas	0.080	11:06
07 Room Air	0.000	11:06

Lot No = 14323080A4
Cyl No = 5
Exp Date = 06/05/2025
County = 64 Oper No. = 133654



Operator Signature
EMILY HENLEY

Remarks:

Form 106-I8000