

## ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59282 (08/2018)

PEWT		Infoxilyzer® Serial Number	
as Lot Number 14323080 A I	Gas Cylinder Number	Gas Expiration Date	

2. Perform an ACA Test (Level 1, Function C).

Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).

A. 0.234 AC

B. 0.23 JAC

C. 0.033 AC

3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.

4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.

5. Send the Following to the Crime Laboratory:

A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).

B. ACA Test Record.

Chemical Test Operator Signature		Date 4-3-24/
Reviewed By (Crime Laboratory Use Only)	Male	Date 08April2024

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = NEWT
04/03/2024

SN 80-006691 8164.14.00 09/16 19:48

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	19:49
02 Std. Gas	0.079	19:49
03 Room Air	0.000	19:50
04 Std. Gas	0.078	19:50
05 Room Air	0.000	19:51
06 Std. Gas	0.079	19:51
07 Room Air	0.000	19:51

Lot No = 14323080A1

Cyl No = 71

Exp Date = 06/05/2025

County = 60

Oper No. = 132000

Operator Signature MARK NEW HOLY

Remarks: ALA Test

Form 106-I8000