

ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

Chemical Test Operator Name (F Kari Stuart	Print)			
Location 1201 12th St SE, Suite C. Watford City, ND 58854		Intoxilyzer® Serial Number 80-003066		
Gas Lot Number 01524080A4	Gas Cylinder Nu	00-003000	00-003066	
71024060A4	17	mber	Gas Expiration Date 03/05/2026	
neck When Done:				

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function 8).
- 2. Perform an ACA Test (Level 1, Function C). Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. O. 18 1 AC B. 0.012 AC C. 0.001 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- ★ 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature		Date 2/5/2025
Reviewed By (Crime Laboratory Use Only) Anna Narehood	Male	05Feb2025

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = MCSO 02/05/2025

Alcohol Analyzer SN 80-003066 8164.14.00 09/16 08:41

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:41
02 Std. Gas	0.081	08:42
03 Room Air	0.000	08:42
04 Std. Gas	0.079	08:42
05 Room Air	0.000	08:43
06 Std. Gas	0.081	08:43
07 Room Air	0.000	08:44

Lot No = 01524080A4

Cyl No = 17

Exp Date = 03/05/2026

County = 27

Oper No. = 132007

Operator Signature KARI STUART

Remarks:

Form 106-I8000

STATE OF NORTH DAKOTA } ...

I certify that the foregoing is a true and correct copy of the original document on file in this office.

Dated at Watford City, North Dayota.