



**ETHANOL GAS STANDARD CYLINDER REPORT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>William Earl</i>		
Location <i>DEVL</i>	Intoxilyzer® Serial Number <i>80-006671</i>	
Gas Lot Number <i>25024080A3</i>	Gas Cylinder Number <i>04</i>	Gas Expiration Date <i>10/05/2026</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).  
 A. 0.*082* AC  
 B. 0.*083* AC  
 C. 0.*082* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.

- Rec'd  
AEN  
27Feb25
- 5. Send the Following to the Crime Laboratory:  
 A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).  
 B. ACA Test Record.

Chemical Test Operator Signature <i>William Earl 9534</i>	Date <i>2/26/2025</i>
Reviewed By (Crime Laboratory Use Only) Anna Narehood	Date 27Feb2025


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006671  
Location = DEVL      8164.14.00 09/16  
02/26/2025      22:53

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	22:54
02 Std. Gas	0.082	22:54
03 Room Air	0.000	22:55
04 Std. Gas	0.083	22:55
05 Room Air	0.000	22:55
06 Std. Gas	0.082	22:56
07 Room Air	0.000	22:56

Lot No = 25024080A3  
Cyl No = 4  
Exp Date = 10/05/2026  
County = 36      Oper No. = 132008



Operator Signature  
WILLIAM EARL

Remarks:

Form 106-I8000