



ETHANOL GAS STANDARD CYLINDER REPORT
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
 SFN 59282 (08/2018)

Chemical Test Operator Name (Print) CHRISTOPHER HUARD # 132570		
Location RICHLAND COUNTY JAZZ	Intoxilyzer® Serial Number 80-003065	
Gas Lot Number Z84Z308DA3	Gas Cylinder Number 47	Gas Expiration Date 11-05-2025

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 A. **0080** AC
 B. **0078** AC
 C. **0078** AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 B. ACA Test Record.

Chemical Test Operator Signature 	Date 02/26/25
Reviewed By (Crime Laboratory Use Only) Anna Narehood	Date 13Mar2025


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-003065
Location = WAHP 8164.14.00 09/16
02/26/2025 13:08

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:09
02 Std. Gas	0.080	13:10
03 Room Air	0.000	13:10
04 Std. Gas	0.078	13:10
05 Room Air	0.000	13:11
06 Std. Gas	0.078	13:11
07 Room Air	0.000	13:12

Lot No = 28423080A3
Cyl No = 47
Exp Date = 11/05/2025
County = 39 Oper No. = 132570



Operator Signature
CHRISTOPHER HUARD

Remarks: GOOD ACA

Form 106-I8000